



PROTECFINA

SAFETY AND PERFORMANCE

HOW TO CHOOSE AND USE YOUR GLOVES

“Food contact” certification “Food safe” symbol



For the “food safe” symbol to figure on glove dispenser boxes, the glove must have undergone successful tests for all types of foodstuffs (acetic acid (3%), ethanol (10%) and refined olive oil). If not, restrictions on use must accompany this symbol clearly. This “food safe” symbol is printed on the back of all our glove dispenser boxes. This means that our gloves have been “food contact” certified by independent bodies in accordance with the applicable European regulations and directives. The test reports are at the disposal of our customers.

Advice for safe use

- Chose the type of gloves based on the risk incurred
- Chose the right size
- Inspect the gloves to check there are no apparent flaws
- Do not wear jewellery on your fingers and avoid long or sharp nails
- Only put gloves on dry, clean hands
- Change gloves frequently, especially if they become damaged
- Remove gloves immediately following contact with a dangerous substance
- Wash your hands before and after wearing the gloves


Size chart

Letter	In inches	GB	France	Palm circumference	Hand length
XS	5/6	Extra-Small	Très Petit	130 mm	151 mm
S	6/7	Small	Petit	152 mm	160 mm
M	7/8	Medium	Moyen	178 mm	171 mm
L	8/9	Large	Grand	203 mm	182 mm
XL	9/10	Extra-Large	Très Grand	229 mm	192 mm

1 inch = 2.54 cm



**Palm
circumference**



**Hand
length**





PROTECFINA
SAFETY AND PERFORMANCE

HOW TO TAKE OFF DISPOSABLE GLOVES



Pinch and hold the outside of the cuff next to your palm



Peel the glove off in one movement, turning it inside out




Keep the inside-out glove rolled up into a ball in the gloved hand



Take hold of the remaining glove with your un-gloved hand



Peel the second glove off over the balled-up glove in your hand



Throw away both gloves in an appropriate container



Wash your hands thoroughly

Recommendations

- Use the right size
- Do not wear jewellery on your fingers
- Avoid false or long nails
- Put gloves on clean, dry hands
- Change gloves frequently
- Avoid excessive stretching